



Beta-lactamase inhibitory protein Antibody, Biotin conjugated

Product Code	CSB-PA327444LD01FOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35804
Immunogen	Recombinant Streptomyces clavuligerus Beta-lactamase inhibitory protein (37-201AA)
Raised In	Rabbit
Species Reactivity	Streptomyces clavuligerus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-lactamase inhibitory protein (BLIP)
Immunogen Species	Streptomyces clavuligerus
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.