



glcB Antibody, Biotin conjugated

Product Code	CSB-PA327661LD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37330
Immunogen	Recombinant Escherichia coli Malate synthase G protein (2-723AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Malate synthase G (MSG) (EC 2.3.3.9), glcB, glc
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	glcB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.