



glcB Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA327661LD01ENV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P37330 |
| Immunogen | Recombinant Escherichia coli Malate synthase G protein (2-723AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Malate synthase G (MSG) (EC 2.3.3.9), glcB, glc |
| Immunogen Species | Escherichia coli (strain K12) |
| Research Area | Others |
| Target Names | glcB |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |