




# rbcl Antibody

<b>Product Code</b>	CSB-PA328678LA01GGV
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P27066
<b>Immunogen</b>	Recombinant Soybean Ribulose biphosphate carboxylase large chain protein (3-475AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Soybean
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ribulose biphosphate carboxylase large chain (RuBisCO large subunit) (EC 4.1.1.39), rbcl
<b>Immunogen Species</b>	Soybean
<b>Research Area</b>	Others
<b>Target Names</b>	rbcl

## Image



### Western Blot

Positive WB detected in Recombinant protein

All lanes: rbcl antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 69 kDa

Observed band size: 69 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.