



# rbcl Antibody, HRP conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA328678LB01GGV  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                            |
| <b>Uniprot No.</b>         | P27066   |
| <b>Immunogen</b>           | Recombinant Soybean Ribulose biphosphate carboxylase large chain protein (3-475AA)       |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Soybean  |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | HRP  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4         |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Ribulose biphosphate carboxylase large chain (RuBisCO large subunit) (EC 4.1.1.39), rbcl |
| <b>Immunogen Species</b>   | Glycine max (Soybean) (Glycine hispida)  |
| <b>Research Area</b>       | Others   |
| <b>Target Names</b>        | rbcl   |