



rbcl Antibody, Biotin conjugated

Product Code	CSB-PA328678LD01GGV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27066
Immunogen	Recombinant Soybean Ribulose biphosphate carboxylase large chain protein (3-475AA)
Raised In	Rabbit
Species Reactivity	Soybean
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribulose biphosphate carboxylase large chain (RuBisCO large subunit) (EC 4.1.1.39), rbcl
Immunogen Species	Glycine max (Soybean) (Glycine hispida)
Research Area	Others
Target Names	rbcl
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.