



prgJ Antibody, HRP conjugated

Product Code	CSB-PA331347LB01SXB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41785
Immunogen	Recombinant Salmonella typhimurium Protein PrgJ protein (1-101AA)
Raised In	Rabbit
Species Reactivity	Salmonella typhimurium
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein PrgJ, prgJ
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Research Area	Others
Target Names	prgJ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.