



# SRP72 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA334457HB01SVG
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P38688
<b>Immunogen</b>	Recombinant Saccharomyces cerevisiae Signal recognition particle subunit SRP72 protein (1-180AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Saccharomyces cerevisiae
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Signal recognition particle subunit SRP72 (Signal recognition particle 72 kDa protein homolog), SRP72
<b>Immunogen Species</b>	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Research Area</b>	Others
<b>Target Names</b>	SRP72
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.