



SRP72 Antibody, Biotin conjugated

Product Code	CSB-PA334457HD01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38688
Immunogen	Recombinant Saccharomyces cerevisiae Signal recognition particle subunit SRP72 protein (1-180AA)
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal recognition particle subunit SRP72 (Signal recognition particle 72 kDa protein homolog), SRP72
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Research Area	Others
Target Names	SRP72
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.