



# WSS1 Antibody

<b>Product Code</b>	CSB-PA334485LA01SVG
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P38838
<b>Immunogen</b>	Recombinant <i>Saccharomyces cerevisiae</i> DNA-dependent metalloprotease WSS1 protein (1-269AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	<i>Saccharomyces cerevisiae</i>
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	DNA-dependent metalloprotease WSS1 (EC 3.4.24.-) (DNA damage response protein WSS1) (SUMO-dependent isopeptidase WSS1) (Weak suppressor of SMT3 protein 1), WSS1
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Research Area</b>	Others
<b>Target Names</b>	WSS1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.