



WSS1 Antibody, Biotin conjugated

Product Code	CSB-PA334485LD01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38838
Immunogen	Recombinant <i>Saccharomyces cerevisiae</i> DNA-dependent metalloprotease WSS1 protein (1-269AA)
Raised In	Rabbit
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-dependent metalloprotease WSS1 (EC 3.4.24.-) (DNA damage response protein WSS1) (SUMO-dependent isopeptidase WSS1) (Weak suppressor of SMT3 protein 1), WSS1
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Research Area	Others
Target Names	WSS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.