



pbpD Antibody, Biotin conjugated

Product Code	CSB-PA334764LD01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40750
Immunogen	Recombinant Bacillus subtilis Penicillin-binding protein 4 protein (213-450AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Penicillin-binding protein 4 (PBP 4) [Includes: Penicillin-insensitive transglycosylase (EC 2.4.1.129) (Peptidoglycan TGase); Penicillin-sensitive transpeptidase (EC 3.4.16.4) (DD-transpeptidase)], pbpD
Immunogen Species	Bacillus subtilis (strain 168)
Research Area	Others
Target Names	pbpD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.