



Tnni3 Antibody, HRP conjugated

Product Code	CSB-PA341785LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48787
Immunogen	Recombinant Mouse Troponin I, cardiac muscle protein (2-211AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Troponin I, cardiac muscle (Cardiac troponin I), Tnni3
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Tnni3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.