



Tnni3 Antibody, FITC conjugated

Product Code	CSB-PA341785LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48787
Immunogen	Recombinant Mouse Troponin I, cardiac muscle protein (2-211AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Troponin I, cardiac muscle (Cardiac troponin I), Tnni3
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Tnni3