




OSM34 Antibody

Product Code	CSB-PA342308LA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50700
Immunogen	Recombinant Arabidopsis thaliana Osmotin-like protein OSM34 protein (23-244AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Osmotin-like protein OSM34, OSM34
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	OSM34

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: OSM34 antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 44 kDa

observed band size: 44 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.