



OSM34 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA342308LC01DOA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P50700 |
| Immunogen | Recombinant Arabidopsis thaliana Osmotin-like protein OSM34 protein (23-244AA) |
| Raised In | Rabbit |
| Species Reactivity | Arabidopsis thaliana |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Osmotin-like protein OSM34, OSM34 |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Research Area | Others |
| Target Names | OSM34 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |