



dlx2a Antibody

Product Code	CSB-PA342691LA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50574
Immunogen	Recombinant Zebrafish Homeobox protein Dlx2a protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Dlx2a (DLX-2) (Distal-less homeobox gene 2a), dlx2a, dlx2
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Neuroscience
Target Names	dlx2a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.