




ftsZ Antibody

Product Code	CSB-PA343512LA01WAP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50907
Immunogen	Recombinant Wolbachia pipientis Cell division protein FtsZ protein (1-315AA)
Raised In	Rabbit
Species Reactivity	Wolbachia pipientis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division protein FtsZ (Fragment), ftsZ
Immunogen Species	Wolbachia pipientis
Research Area	Others
Target Names	ftsZ

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ftsZ antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 39 kDa

observed band size: 39 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.