



ftsZ Antibody, HRP conjugated

Product Code	CSB-PA343512LB01WAP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50907
Immunogen	Recombinant Wolbachia pipientis Cell division protein FtsZ protein (1-315AA)
Raised In	Rabbit
Species Reactivity	Wolbachia pipientis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division protein FtsZ (Fragment), ftsZ
Immunogen Species	Wolbachia pipientis
Research Area	Others
Target Names	ftsZ