



ftsZ Antibody, Biotin conjugated

Product Code	CSB-PA343512LD01WAP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50907
Immunogen	Recombinant Wolbachia pipientis Cell division protein FtsZ protein (1-315AA)
Raised In	Rabbit
Species Reactivity	Wolbachia pipientis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division protein FtsZ (Fragment), ftsZ
Immunogen Species	Wolbachia pipientis
Research Area	Others
Target Names	ftsZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.