



Probable terminase, ATPase subunit Antibody, FITC conjugated

Product Code	CSB-PA343952LC01HAN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51718
Immunogen	Recombinant Haemophilus phage HP1 Probable terminase, ATPase subunit protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Haemophilus phage HP1
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable terminase, ATPase subunit (ORF16)
Immunogen Species	Haemophilus phage HP1 (strain HP1c1) (Bacteriophage HP1)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.