




petE Antibody

Product Code	CSB-PA344627HA01PSD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50057
Immunogen	Recombinant Prochlorothrix hollandica Plastocyanin protein (35-131AA)
Raised In	Rabbit
Species Reactivity	Prochlorothrix hollandica
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plastocyanin, petE
Immunogen Species	Prochlorothrix hollandica
Research Area	Others
Target Names	petE

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: petE antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 27 kDa

observed band size: 27 kDa