



petE Antibody, Biotin conjugated

Product Code	CSB-PA344627HD01PSD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50057
Immunogen	Recombinant Prochlorothrix hollandica Plastocyanin protein (35-131AA)
Raised In	Rabbit
Species Reactivity	Prochlorothrix hollandica
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plastocyanin, petE
Immunogen Species	Prochlorothrix hollandica
Research Area	Others
Target Names	petE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.