



ponA Antibody, Biotin conjugated

Product Code	CSB-PA345793LD01DKK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54660
Immunogen	Recombinant Slime mold Ponticulin protein (23-118AA)
Raised In	Rabbit
Species Reactivity	Slime mold
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ponticulin, ponA
Immunogen Species	Dictyostelium discoideum (Slime mold)
Research Area	Others
Target Names	ponA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.