




hfb1 Antibody

Product Code	CSB-PA347143HA01HYE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52754
Immunogen	Recombinant Hypocrea jecorina Hydrophobin-1 protein (23-97AA)
Raised In	Rabbit
Species Reactivity	Hypocrea jecorina
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydrophobin-1 (Hydrophobin I) (HFBI), hfb1
Immunogen Species	Hypocrea jecorina
Research Area	Others
Target Names	hfb1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: hfb1 antibody at 4μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 25 kDa

observed band size: 25 kDa