



X Antibody, Biotin conjugated

Product Code	CSB-PA355969LD01HEM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03165
Immunogen	Recombinant Hepatitis B virus genotype D subtype ayw Protein X protein (1-154AA)
Raised In	Rabbit
Species Reactivity	Hepatitis B virus genotype D subtype ayw
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein X (HBx) (Peptide X) (pX), X
Immunogen Species	Hepatitis B virus genotype D subtype ayw (isolate France/Tiollais/1979) (HBV-D)
Research Area	Others
Target Names	X
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.