



# htpG Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA358897LC01ENV
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0A6Z3
<b>Immunogen</b>	Recombinant Escherichia coli Chaperone protein HtpG protein (1-624AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Escherichia coli
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Chaperone protein HtpG (Heat shock protein C62.5) (Heat shock protein HtpG) (High temperature protein G), htpG
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Research Area</b>	Others
<b>Target Names</b>	htpG