



htpG Antibody, Biotin conjugated

Product Code	CSB-PA358897LD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A6Z3
Immunogen	Recombinant Escherichia coli Chaperone protein HtpG protein (1-624AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chaperone protein HtpG (Heat shock protein C62.5) (Heat shock protein HtpG) (High temperature protein G), htpG
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	htpG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.