



C3 Antibody

Product Code	CSB-PA360475HA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01026
Immunogen	Recombinant Rat Complement C3 protein (671-746AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C3 [Cleaved into: Complement C3 beta chain; C3-beta-c (C3bc) (Neutrophil chemotactic factor-2) (ENCF-2); Complement C3 alpha chain; C3a anaphylatoxin (Neutrophil chemotactic factor-1) (ENCF-1); Acylation stimulating protein (ASP) (C3adesArg); Complement C3b alpha' chain; Complement C3c alpha' chain fragment 1; Complement C3dg fragment; Complement C3g fragment; Complement C3d fragment; Complement C3f fragment; Complement C3c alpha' chain fragment 2], C3
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Immunology
Target Names	C3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.