



# C3 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA360475HB01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P01026
<b>Immunogen</b>	Recombinant Rat Complement C3 protein (671-746AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Complement C3 [Cleaved into: Complement C3 beta chain; C3-beta-c (C3bc) (Neutrophil chemotactic factor-2) (ENCF-2); Complement C3 alpha chain; C3a anaphylatoxin (Neutrophil chemotactic factor-1) (ENCF-1); Acylation stimulating protein (ASP) (C3adesArg); Complement C3b alpha' chain; Complement C3c alpha' chain fragment 1; Complement C3dg fragment; Complement C3g fragment; Complement C3d fragment; Complement C3f fragment; Complement C3c alpha' chain fragment 2], C3
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Immunology
<b>Target Names</b>	C3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.