




Lysozyme C-1 Antibody

Product Code	CSB-PA360572LA01BZD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00705
Immunogen	Recombinant Anas platyrhynchos Lysozyme C-1 protein (19-147AA)
Raised In	Rabbit
Species Reactivity	Anas platyrhynchos
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysozyme C-1 (EC 3.2.1.17) (1,4-beta-N-acetylmuramidase C)
Immunogen Species	Anas platyrhynchos
Research Area	Others

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: Lysozyme C-1 antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.