



Cationic trypsin Antibody, HRP conjugated

Product Code	CSB-PA361971LB01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06871
Immunogen	Recombinant Dog Cationic trypsin protein (24-246AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cationic trypsin (EC 3.4.21.4)
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.