



# Cationic trypsin Antibody, HRP conjugated

Product Code	CSB-PA361971LB01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06871
Immunogen	Recombinant Dog Cationic trypsin protein (24-246AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cationic trypsin (EC 3.4.21.4)
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Research Area	Others