



mrcA Antibody, HRP conjugated

Product Code	CSB-PA363014LB01NGG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A0Z6
Immunogen	Recombinant Neisseria meningitidis serogroup B Penicillin-binding protein 1A protein (206-413AA)
Raised In	Rabbit
Species Reactivity	Neisseria meningitidis serogroup B
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Penicillin-binding protein 1A (PBP-1a) (PBP1a) [Includes: Penicillin-insensitive transglycosylase (EC 2.4.1.129) (Peptidoglycan TGase); Penicillin-sensitive transpeptidase (EC 3.4.16.4) (DD-transpeptidase)], mrcA, ponA
Immunogen Species	Neisseria meningitidis serogroup B (strain MC58)
Research Area	Others
Target Names	mrcA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.