



# mrcA Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA363014LD01NGG
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0A0Z6
<b>Immunogen</b>	Recombinant Neisseria meningitidis serogroup B Penicillin-binding protein 1A protein (206-413AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Neisseria meningitidis serogroup B
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Penicillin-binding protein 1A (PBP-1a) (PBP1a) [Includes: Penicillin-insensitive transglycosylase (EC 2.4.1.129) (Peptidoglycan TGase); Penicillin-sensitive transpeptidase (EC 3.4.16.4) (DD-transpeptidase)], mrcA, ponA
<b>Immunogen Species</b>	Neisseria meningitidis serogroup B (strain MC58)
<b>Research Area</b>	Others
<b>Target Names</b>	mrcA
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.