



# mscL Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA363948LB01EOD
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0A743
<b>Immunogen</b>	Recombinant Escherichia coli O157:H7 Large-conductance mechanosensitive channel protein (1-136AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Escherichia coli O157:H7
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Large-conductance mechanosensitive channel, mscL
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Research Area</b>	Others
<b>Target Names</b>	mscL
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.