



mscL Antibody, Biotin conjugated

Product Code	CSB-PA363948LD01EOD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A743
Immunogen	Recombinant Escherichia coli O157:H7 Large-conductance mechanosensitive channel protein (1-136AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli O157:H7
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Large-conductance mechanosensitive channel, mscL
Immunogen Species	Escherichia coli O157:H7
Research Area	Others
Target Names	mscL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.