



# clpP Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA401627LB01CWV
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	A5I6W1
<b>Immunogen</b>	Recombinant Clostridium botulinum ATP-dependent Clp protease proteolytic subunit protein (1-194AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Clostridium botulinum
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	ATP-dependent Clp protease proteolytic subunit (EC 3.4.21.92) (Endopeptidase Clp), clpP
<b>Immunogen Species</b>	Clostridium botulinum (strain Hall / ATCC 3502 / NCTC 13319 / Type A)
<b>Research Area</b>	Others
<b>Target Names</b>	clpP
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.