



clpP Antibody, Biotin conjugated

Product Code	CSB-PA401627LD01CWV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5I6W1
Immunogen	Recombinant Clostridium botulinum ATP-dependent Clp protease proteolytic subunit protein (1-194AA)
Raised In	Rabbit
Species Reactivity	Clostridium botulinum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent Clp protease proteolytic subunit (EC 3.4.21.92) (Endopeptidase Clp), clpP
Immunogen Species	Clostridium botulinum (strain Hall / ATCC 3502 / NCTC 13319 / Type A)
Research Area	Others
Target Names	clpP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.