



porB Antibody, FITC conjugated

Product Code	CSB-PA518750LC01NGH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E6MZM0
Immunogen	Recombinant Neisseria meningitidis serogroup B / serotype 15 Major outer membrane protein P.IB protein (20-331AA)
Raised In	Rabbit
Species Reactivity	Neisseria meningitidis serogroup B / serotype 15
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major outer membrane protein P.IB (PIB) (Protein IB) (Class 3 protein) (Porin), porB
Immunogen Species	Neisseria meningitidis serogroup B / serotype 15 (strain H44/76)
Research Area	Others
Target Names	porB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.