



# porB Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA518750LD01NGH
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	E6MZM0
<b>Immunogen</b>	Recombinant Neisseria meningitidis serogroup B / serotype 15 Major outer membrane protein P.IB protein (20-331AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Neisseria meningitidis serogroup B / serotype 15
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Major outer membrane protein P.IB (PIB) (Protein IB) (Class 3 protein) (Porin), porB
<b>Immunogen Species</b>	Neisseria meningitidis serogroup B / serotype 15 (strain H44/76)
<b>Research Area</b>	Others
<b>Target Names</b>	porB
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.