



yvgO Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA520104LC01BRJ |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O32211 |
| Immunogen | Recombinant Bacillus subtilis Stress response protein YvgO protein (27-161AA) |
| Raised In | Rabbit |
| Species Reactivity | Bacillus subtilis |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Stress response protein YvgO, yvgO |
| Immunogen Species | Bacillus subtilis (strain 168) |
| Research Area | Others |
| Target Names | yvgO |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |