



yvgO Antibody, Biotin conjugated

Product Code	CSB-PA520104LD01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O32211
Immunogen	Recombinant Bacillus subtilis Stress response protein YvgO protein (27-161AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stress response protein YvgO, yvgO
Immunogen Species	Bacillus subtilis (strain 168)
Research Area	Others
Target Names	yvgO