



yvgO Antibody, Biotin conjugated

Product Code	CSB-PA520104LD01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O32211
Immunogen	Recombinant Bacillus subtilis Stress response protein YvgO protein (27-161AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stress response protein YvgO, yvgO
Immunogen Species	Bacillus subtilis (strain 168)
Research Area	Others
Target Names	yvgO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.