




Fgf15 Antibody

Product Code	CSB-PA522052LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35622
Immunogen	Recombinant Mouse Fibroblast growth factor 15 protein (26-218AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Fgf15

Image



Western Blot

Positive WB detected in: Mouse Liver tissue lysate, Mouse Spleen tissue lysate, Mouse Lung tissue

All lanes: Fgf15 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 26 kDa

Observed band size: 26 kDa