




SAP9 Antibody

Product Code	CSB-PA522117LA01CZF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O42779
Immunogen	Recombinant Yeast Candidapepsin-9 protein (18-544AA)
Raised In	Rabbit
Species Reactivity	Yeast
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Candidapepsin-9 (EC 3.4.23.24) (ACP 9) (Aspartate protease 9) (Secreted aspartic protease 9), SAP9
Immunogen Species	Yeast
Research Area	Others
Target Names	SAP9

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: SAP9 antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 75 kDa

observed band size: 75 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.