



# SAP9 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA522117LB01CZF
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O42779
<b>Immunogen</b>	Recombinant Yeast Candidapepsin-9 protein (18-544AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Yeast
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Candidapepsin-9 (EC 3.4.23.24) (ACP 9) (Aspartate protease 9) (Secreted aspartic protease 9), SAP9
<b>Immunogen Species</b>	Candida albicans (Yeast)
<b>Research Area</b>	Others
<b>Target Names</b>	SAP9
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.