



vegfaa Antibody, Biotin conjugated

Product Code	CSB-PA524826LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O73682
Immunogen	Recombinant Zebrafish Vascular endothelial growth factor A-A protein (24-188AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vascular endothelial growth factor A-A (VEGF-A-A), vegfaa, vegf vegfa
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	vegfaa
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.