



ELP6 Antibody, HRP conjugated

Product Code	CSB-PA608122LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0PNE2
Immunogen	Recombinant Human Elongator complex protein 6 protein (95-195AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongator complex protein 6 (Angiotonin-transactivated protein 1) (Protein TMEM103), ELP6, ATP1 C3orf75 TMEM103
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ELP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.