



ELP6 Antibody, HRP conjugated

Product Code	CSB-PA608122LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0PNE2
Immunogen	Recombinant Human Elongator complex protein 6 protein (95-195AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongator complex protein 6 (Angiotonin-transactivated protein 1) (Protein TMEM103), ELP6, ATP1 C3orf75 TMEM103
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ELP6